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Inflation Device Instruction for Use

【Precaution】

- ▶ Only operate by professional physician, please read the IFU carefully before operating the device.
- ▶ Do not use, if the package is open, break, leak or expired the effective period of sterilization.
- ▶ This device was sterilized by EO and disposable, RSTERILIZATION was prohibited.
- ▶ Refer to the Limited Burst Pressure of balloon dilatation catheter or other device connecting to the device provided by original manufacture, before operating.
- ▶ Loss of pressure may indicate a leak in the system.


【Description】

The Inflation Device is a 20mL disposable device capable of producing a maximum pressure of 30 ATM and 40 ATM. The device is composed of handle, top and bottom shells, push block, syringe, pressure gage, pressure tubing and male luer locker.

【Variants】

Variant	Explanation
K-P-30	Maximum Pressure 30ATM, Normal handle, Without three-way connection
K-P-30S	Maximum Pressure 30ATM, Normal handle, With three-way connection
K-P-30I	Maximum Pressure 30ATM, I shape handle, Without three-way connection
K-P-30IS	Maximum Pressure 30ATM, I shape handle, With three-way connection
K-P-40	Maximum Pressure 40ATM, Normal handle, Without three-way connection
K-P-40S	Maximum Pressure 40ATM, Normal handle, With three-way connection
K-P-40I	Maximum Pressure 40ATM, I shape handle, Without three-way connection
K-P-40IS	Maximum Pressure 40ATM, I shape handle, With three-way connection

Note: K: abbreviation of Kossel; P: refer to Inflation Device; I: I shape handle; S: with a three-way connection

【Intended Use】

The inflation device is used to inflate or deflate an angioplasty balloon, and to measure the pressure within the balloon.

【Instruction for Use】

- Note: Before operating, assure the package and production is unbroken.
- Let both of the push block and pull rod backward at the same time. Release the push block when the volume of dilution contrast or other liquid absorbing by the device meet the suitable quantity.
 - Keep the device upward, let the push block backward and let the pull rod forward, exclude the residual air in the syringe. If necessary, release the push block and rotate the pull rod suitably, exclude the residual air in the syringe.
 - Connect the luer locker to the balloon catheter, keep the graduation line facing to the operator.
 - Let both of the push block and pull rod backward, exclude the residual air in the syringe.
 - Repeat the second step to exclude the residual air in the syringe.

● Release the push block, inflate the balloon catheter by rotating the pull rod clockwise, and deflate the balloon catheter by rotating the pull rod anti-clockwise. When the pressure decline to no high than 10atm, the pressure of the balloon catheter can be released by making the push block backward.

● Let both of the push block and pull rod backward, and release the push block, the balloon catheter will be vacuum.

● The pressure could loss quickly by making the push block backward, when the inflation device is keeping pressure.

Note: Risks may occur when release the pressure too quickly, the operator should operate suitably according to clinical situation.

【Storage】
The product should be stored in a environment where the relative humidity is not more than 80%, the temperature is below 40°C, and there is no corrosive gas, cool, dry, well ventilated and clean environment.

【Validity】
Two years

【Explanation of Symbols】

- | | |
|--|--|
|  Do not reuse |  Use by |
|  Catalog Number |  Manufacture Date |
|  Lot number |  Sterilized using Ethylene Oxide |
|  Consult instructions for use |  Do not use if package is damaged |
|  Keep Away from Rain |  Keep Away from Sunlight |
|  This Way Up |  Caution |
|  Temperature Limits |  Fragile |

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Kossel Medtech (Suzhou) Co., Ltd.
Name: MedPath GmbH